

Identification and characterization of a QTL for growth of *Fusarium circinatum* on pine-based medium

Benedicta S. Swalarsk-Parry¹, E.T. Steenkamp¹, S. van Wyk¹, Q.C. Santana^{1,2}, M.A. van der Nest^{1,2}, A. Hammerbacher³, B.D. Wingfield¹ and L. De Vos¹

¹Department of Biochemistry, Genetics and Microbiology, Forestry and Agricultural Biotechnology Institute (FABI), University of Pretoria, Pretoria 0002, South Africa.

²Agricultural Research Council, Biotechnology Platform, 100 Old Soutpan Road, Onderstepoort, Pretoria, South Africa.

³Department of Zoology and Entomology, Forestry and Agricultural Biotechnology Institute (FABI), University of Pretoria, Pretoria 0002, South Africa.

Supplementary File S2

The sequence of the primers used to confirm the absence or presence of genes and regions targeted in this study.

Primer name	Primer sequence (5' to 3')	Target gene	Target region
P _{73F}	ACTCCCCAACTGACTCATGG	g13407	1 194 142 - 1 195 143
P _{73R}	CCACTCGCTTCCAAACACAA		
P _{74F}	CAGTTGGCACTAGTCCTCCA	g13408	1 198 252 - 1 198 921
P _{74R}	TCAGGATCGATACCAAGGGC		
P _{4F}	TTGCCCTCTCCAGTACCAAG	FSP34 indel	1 119 532 - 1 127 164
P _{5R}	TGTCTTTGGTTCACACGCTG		
KS17 _F	GATGACAAACCGGCTCACTC	g13276	1 206 036 - 1 208 265
KS17 _R	CAGCCATGAGACGGTAAAGC		
FSP34 _F	CATGGACGCGTCAAGTACTG	g13376	1 165 752 - 1 168 139
FSP34 _R	TGCTTGCAGGTTCAAAGGTC		
FT _F	GGTGGTGCTAAGGACTCACT	g13420	1 238 853 - 1 242 201
FT _R	GGTGCTCGCTTTCCAATAGG		

Cycling conditions using the KAPA Taq PCR kit (Kapa Biosystems) – these conditions were mainly used for the genes g13407 and g13408

Step	Temperature	Time	Cycles
Initial Denaturation	95 °C	3 min	1
Denaturation	95 °C	30 seconds	30
Annealing	variable *	30 seconds	
Extension	72 °C	1 min per kb	1
Final Extension	72 °C	5 minutes	
Hold	4 °C	∞	

*Annealing temperature for primer set P_{73F+R} = 57 °C primer set P_{74F+R} = 56 °C

Cycling conditions using the LongAmp® Taq 2X Master Mix (New England BioLabs) – these conditions were mainly used for the indel in *F. circinatum* FSP34

Step	Temperature	Time	Cycles
Initial Denaturation	94 °C	30 seconds	1
Denaturation	94 °C	30 seconds	
Annealing	55 °C	1 minute	10
Extension	65 °C	8 minutes	
Denaturation	94 °C	30 seconds	
Annealing	55 °C	1 minute	20
Extension	65 °C	8 minutes + 20 sec each cycle	
Final Extension	65 °C	10 minutes	1
Hold	4 °C	∞	

LongAmp® Taq 2X Master Mix (New England BioLabs) protocol to amplify three different genes* in CMWF350 (FSP34), CMWF389, and CMWF674 (KS17).

Step	Temperature	Time	Cycles
Initial Denaturation	94 °C	30 seconds	1
Denaturation	94 °C	30 seconds	
Annealing	52 °C	30 seconds	30
Extension	65 °C	2 minutes 30 seconds	
Final Extension	65 °C	10 minutes	1
Hold	4 °C	∞	

*Primer set FSP34_{F+R}, KS17_{F+R}, and FT_{F+R} were used to amplify the different genes in *F. circinatum* (FSP34 and KS17) and *F. temperatum* (CMWF389).