

# ANTIMICROBIAL EFFECTIVENESS OF ELECTRO-CHEMICALLY ACTIVATED WATER AS AN ENDODONTIC IRRIGATION SOLUTION

by

WAYNE PHILLIP WILLIAMS

BDS (Witwatersrand)

Supervisor: DR JT MARAIS

Submitted in partial fulfillment of the requirements for the degree

MChD (Prosthodontics)

in the

FACULTY OF DENTISTRY
UNIVERSITY OF PRETORIA
PRETORIA

June 1999



# **DECLARATION**

I hereby declare that this dissertation, submitted by me in partial fulfillment of the requirements for the degree of MChD (Prosthodontics) at the University of Pretoria, has not previously been submitted for a degree at any other University

WAYNE P WILLIAMS

DATE



# **Declaration by Supervisor**

As the supervisor of this dissertation submitted by **Dr WP Williams** entitled:

# ANTIMICROBIAL EFFECTIVENESS OF ELECTRO-CHEMICALLY ACTIVATED WATER AS AN ENDODONTIC IRRIGATION SOLUTION

I approve the document for submission and further declare that an article for publication in the *International Endodontic Journal* has been prepared

This dissertation has not previously been submitted for any degree at any other university

This dissertation complies with all the rules as set out in the University of Pretoria yearbook (Reg G24(d); G24(e); G25; G32(e); G33; G37)

Dr JT Marais

Date



### **ACKNOWLEDGEMENTS**

My sincere gratitude is expressed to **Dr JT Marais** for his continued enthusiasm, encouragement and intense devotion to the cause of this research. His integrity and irreproachable manner of conduct throughout this project has served as a fine example

I am indebted to the late **Prof WJC Coetzee** for his scientific accuracy and unselfish dedication to the correctness of this manuscript

I wish to thank the Centre for Stomatological Research – University of Pretoria, in particular, Dr Francien Botha for assistance with microbiological techniques as well as Ms I De Wet and Ms A Grobler for their unselfish devotion of time beyond their normal line of duty

Through his infinite sense of diplomacy and wisdom, **Prof LH Becker** remains a constant influence on my development as a young and inexperienced clinician and researcher in the field of prosthodontics. The belief and trust that he has placed in me, allowing me the freedom to experiment and experience all that has been in the past forty months, is certainly encouragement for any student from his mentor

**Prof CH Driessen**, an experienced and thorough researcher has been a constant source of positive criticism and encouragement. His assistance with SEM analysis is also appreciated

**Prof PL Kemp**, thank you for the way you have taught me to "think" about research and reporting on an international platform. "There are no answers, only interpretations"

Ms Laura Matzner did statistical revision. I am most grateful for her professional and prompt service

I will always be indebted to my family (**Henriëtte**, **Warren** and **Rowan**) for allowing me the time away from them to complete this work



### **SUMMARY**

The use of sodium hypochlorite (NaOCl) as an endodontic irrigation solution is effective in eliminating microorganisms from the root canals of human teeth. The purpose of this in vitro study was to evaluate the antimicrobial effectiveness of electro-chemically activated water (ECA) as an endodontic irrigation solution. ECA is a solution comprising many different ions, molecules and radicals. Sixty caries-free, single root, adult, maxillary, anterior human teeth were used. The root canals were instrumented and irrigated with NaOCl in a similar method to that employed for in vivo root canal treatment. The external root surface of each tooth was sealed and the access cavity kept patent so that root canals could be inoculated with a suspension containing four bacteria. The teeth were randomly divided into four groups (n=15). Each group was irrigated ultrasonically, using solutions of distilled water (control), NaOCl (3.5%), and ECA, the latter at pH's 7.0 and 9.0 respectively. Antimicrobial effectiveness was established directly after irrigation and again seven days later, by counting colony forming units on blood agar plates and by spectrophotometric analysis. The surfaces of the root canals were visually examined by means of scanning electron microscopy. Large numbers of bacteria were present in the canals of teeth irrigated with distilled water. No bacteria were observed following irrigation with NaOCl. Neither of the ECA solutions were found to be effective against all the bacteria. Although some reduction in the number of bacteria was evident in the ECA groups, this was not statistically significant (p>0.05). Within the confines of this study ECA did not demonstrate antimicrobial effectiveness.

<sup>&</sup>lt;sup>1</sup> Prevotella intermedia
Porphyromonas gingivalis
Enterococcus faecalis
Actinobacillus actinomycetemcomitans



### **OPSOMMING**

Die gebruik van natriumhipochloriet (NaOCl) as 'n endodontiese irrigeringsmiddel elimineer mikrorganismes in die wortelkanaalstelsel van menslike tande doeltreffend. Die doel van hierdie in vitro studie was om die antimikrobiale effektiwiteit van elektrochemies geaktiveerde water (ECA) as 'n endodontiese irrigasie middel te evalueer. ECA is 'n oplossing van verskeie ione, molekules en radikale. Sestig karies-vrye, enkelkanaal, volwasse, maksillêre, anterior menslike tande is gebruik in hierdie studie. Die wortelkanale is voorberei en geirrigeer met NaOCl volgens die metode wat gebruik word vir in vivo wortelkanaal behandeling. Die eksterne wortel oppervlak van elke tand is geseël en die toegangsopening oopgehou sodat die wortelkanale geinokuleer kon word met 'n suspensie van vier bakterië<sup>2</sup>. Die tande is lukraak in vier groepe (n=15) verdeel. Die vier groepe is elk ultrasonies geirrigeer – die eerste groep met gedistilleerde water (kontrole); die tweede met NaOCl (3.5%) en die derde en vierde met ECA van pH 7.0 en 9.0 onderskeidelik. Antimikrobiale effektiwiteit is vasgestel direk na irrigasie en weer na sewe dae, deur tellings van kolonievormende eenhede op bloed agar plate en deur spektrofotometriese analise. Die dentinale oppervlaktes van die wortelkanaal is ondersoek met die gebruik van skandeer elektron mikroskopie. Groot getalle bakterië was teenwoordig in die kanale van die tande wat met gedistilleerde water geïrrigeer is. Geen bakterië kon waargeneem word na irrigasie met NaOCl nie. Nie een van die ECA oplossings was effektief teen al die bakterië nie. Alhoewel 'n vermindering in die getal bakterië in die ECA groepe gevind is, was dit nie statisties betekenisvol nie (p>0.05). Binne die beperkinge van hierdie studie was daar geen duidelike bewys van die antimikrobiale effektiwiteit van ECA nie.

<sup>&</sup>lt;sup>2</sup> Prevotella intermedia Porphyromonas gingivalis Enterococcus faecalis Actinobacillus actinomycetemcomitans



# LIST OF CONTENTS

CHAPTER I	
Introduction	1
CHAPTER 2	
Literature Review	2
2.1 Technical and biological problems with the use of sodium	
hypochlorite	2
2.2 Alternatives for sodium hypochlorite	3
CHAPTER 3	
Purpose	8
CHAPTER 4	
Materials and Methods	8
4.1 Preparation of teeth	9
4.2 Preparation for first microbiological analysis	10
4.3 Irrigation of teeth	14
4.4 Second microbiological analysis	14
4.5 Statistical analysis	15
4.6 Scanning electron microscopy	15
CHAPTER 5	
Results	17



# CHAPTER 6 Discussion 26 CHAPTER 7 Conclusion 29 REFERENCES 30 ADDENDA 34



# LIST OF FIGURES

Figu	ıre	
1	SEM - Irrigation with NaOCl (x500)	3
2	SEM - Irrigation with ECA (x500)	4
3	SEM - Irrigation with NaOCl (x2500)	4
4	SEM - Irrigation with ECA (x2500)	5
5	STEDS unit for the production of ECA	6
6	Blood agar plate from pilot study to demonstrate antimicrobial	
	properties of ECA	7
7	Higher magnification of blood agar plate seen in Fig.6	7
8	Teeth used for this study	8
9	Spiral filler introduction of calcium hydroxide (Ca(OH)2)	
	(note the sealed external root surface)	10
10	Sterile pipette for obtaining samples from root canal	11
11	Placing sampled solution onto the blood agar plate	11
12	Standardized suspensions of bacteria (Mac Farland-1)	12
13	Anaerobic bacteria imported from ATCC (USA)	12
14	Tooth removed from BHI broth and held upright	
	(note BHI broth in beaker)	13
15	SEM -no irrigation (x6000)	21
16	SEM - distilled water (x6000)	22
17	SEM - NaOCl (x1500)	22
18	SEM - NaOCl (x6000)	23
19	SEM - ECA pH 7.0 (x1500)	23
20	SEM - ECA pH 7.0 (x6000)	24
21	SEM - ECA pH 9.0 (x1500)	24
22	SEM – ECA pH 9.0 (x6000)	25



# LIST OF TABLES

Irrigation solutions and Groups	14
Stages at which teeth were removed from the study for	
SEM examination	16
Colony forming units and spectrophotometric values measured	
immediately after irrigation of root canals with four different	
solutions	17
Colony forming units and spectrophotometric values measured	
seven days after irrigation of root canals with four different	
solutions	18
Increase in colony forming units and spectrophotometric values	
for the period directly following irrigation to the seventh day	
post-irrigation (Table 4 minus Table 3)	19
Dunnett's test results on spectrophotometric values, comparing	
each of Groups B, C and D to Group A (control)	20
Results of paired t-test on spectrophotometric values for bacterial	
growth over the seven day period after irrigation with four	
different solutions	20
	Stages at which teeth were removed from the study for SEM examination  Colony forming units and spectrophotometric values measured immediately after irrigation of root canals with four different solutions  Colony forming units and spectrophotometric values measured seven days after irrigation of root canals with four different solutions  Increase in colony forming units and spectrophotometric values for the period directly following irrigation to the seventh day post-irrigation (Table 4 minus Table 3)  Dunnett's test results on spectrophotometric values, comparing each of Groups B, C and D to Group A (control)  Results of paired t-test on spectrophotometric values for bacterial growth over the seven day period after irrigation with four



# LIST OF ADDENDA

Appendix A	Colony forming units and spectrophotometric analysis	š
	measured directly after irrigation	34
Appendix B	Colony forming units and spectrophotometric analysis	;
	measured seven days after irrigation	35
Appendix C	Molecules and ions contained in ECA	36



### **CHAPTER 1**

### Introduction

The treatment of choice for periapical disease is elimination of microorganisms and their products from the root canal system by means of a root canal treatment. The ADA (American Dental Association) estimates that by the year 2000, thirty million root canals will be treated annually in the USA alone (American Association of Endodontics recertification document, (1989)).

The procedure of root canal treatment in human teeth is often described by the following three terms: Cleaning, Shaping, Filling. The "cleaning" phase consists essentially of the irrigation of the canal with a liquid agent (irrigation solution). The objectives of this phase are destruction of microorganisms; dissolution of vital and necrotic tissue; flushing away of organic debris and bacterial remnants; lubrication of endodontic instruments; prevention of blockage of canals by the by-products of the "shaping" phase, and removal of the smear layer. The focus of this study is on the destruction of microorganisms within the root canal system.

NaOCl at a concentration of between 2.5 and 5.25% is universally accepted as the irrigant of choice for root canal therapy (Cheung and Stock, 1993; Gulabivala and Stock, 1995; Walton and Rivera, 1996).

NaOCl is delivered into the root canal by means of a syringe and needle or an ultrasonic unit. The use of an ultrasonic unit as an irrigation instrument is considered to be superior to the conventional syringe and needle (Cheung and Stock, 1993; Cameron, 1995).

NaOCl achieves the above objectives of the "cleaning" phase with remarkable efficiency, but is also associated with and may cause many complications and problems, some of which may even be life threatening (Jeansonne and White, 1994).



### **CHAPTER 2**

### **Literature Review**

### 2.1 Technical and biological problems with the use of sodium hypochlorite

NaOCl is extremely toxic to living tissue. Extrusion of NaOCl through the apices of teeth causes severe pain, swelling and necrosis (Brown *et al.*, 1995). Accidental extrusion of NaOCl into the maxillary sinus may cause serious complications and even death (Cymbler and Ardakani, 1994). The use of local anaesthetic cartridges filled with NaOCl can also cause severe complications if inadvertently injected into a patient's soft tissue (Caliskan, Turkun and Alper, 1994). Not only is the use of such a cartridge system dangerous, but is also more time consuming as each cartridge must be sterilized and filled.

Because of the corrosive nature of NaOCl, ultrasonic units used for canal irrigation are prone to mechanical blockage and breakdown. Maintenance of such systems is time consuming and expensive. Due to the vile taste and odor of NaOCl, accidental spillage into the mouth during root canal treatment results in severe discomfort for the patient. This invariably requires the patient to rinse the mouth, an action that defeats the object of "cleaning" the root canal, because saliva containing great numbers of bacteria may enter the tooth and root canal during the rinsing process (Yesilsoy *et al.*, 1995).


Because of the above-mentioned problems, dentists often neglect the "cleaning" phase of root canal treatment. This may result in failure of the root canal treatment and a host of other complications such as dental abscesses, root end surgical operations, general anesthesia, hospitalization, further endodontic failures and loss of teeth that had previously been treated.

Substrate adherence capacity of inflammatory macrophages has been shown to decrease significantly with the use of NaOCl (Jimenez-Rubio, *et al*, 1997). Given that adhesion is the first step in the phagocytic process of macrophages and in antigen presentation, NaOCl could inhibit macrophage function and reduce inflammatory reactions in periapical tissues when used in root canal therapy.

### 2.2 Alternatives for sodium hypochlorite

Water, disinfectants and various detergents have been evaluated for possible use as root canal irrigants. Antimicrobial properties of calcium hydroxide (Ca (OH)<sub>2</sub>), chlorhexidine and camphorated paramonochlorophenol were examined by Barbosa *et al.*, (1997). Nd-YAG laser irradiation (Fegan and Steiman, 1995), citric acid (Georgopoulou *et al.*, 1994), oxidative potential water (Hata *et al.*, 1996), chlorhexidine gluconate (Jeansonne and White, 1994) and gluteraldehyde (Jimenez-Rubio *et al.*, 1997) have also been reported on as possible intracanal agents. None of these have been accepted as a replacement for NaOCl for root canal irrigation.

Electro-chemically Activated Water: Russian scientists developed a process whereby they could produce ECA. To date, only one reference to this solution is to be found in western scientific literature (Marais, 1998). Marais irrigated one group of root canals with ECA (catholyte) and another with NaOCl (3.5%), with solutions being delivered from an ultrasonic unit. Root canal surfaces were examined using a scanning electron microscope (Fig 1, 2, 3 & 4). ECA produced surfaces that appeared as clean as those for the NaOCl group. Marais also showed that ECA removed the smear layer, from many surfaces of the canal.








Fig 2: Irrigation with ECA (x500)




Fig 3: Irrigation with NaOCl(x2500)







Fig 4: Irrigation with ECA(x2500)



ECA is the subject of more than 300 Russian patents. More than 20 000 units producing ECA are in operation in Russian hospitals today. The physical and chemical nature of ECA is not yet fully understood. The solution supposedly exists in a metastable or disequilibrious state after production and contains many free radicals and a variety of molecules and ions (see Appendix III, pg. 30). In the metastable state the solution has a very high oxidation-reduction potential. Two types of ECA solution are produced:

- Anolyte, with a high oxidation potential (plus 400 to plus 1200 millivolts). It is
  possible to produce acidic, neutral or alkaline anolyte (pH 2-9). Anolyte is
  claimed to be antimicrobial (according to the manufacturers)<sup>1</sup>.
- Catholyte, an alkaline solution (pH 7-12) with a high reduction potential (minus 80 to minus 900 millivolts, measuring current in opposite direction). Catholyte is reputed to have a strong cleaning or detergent effect.

<sup>&</sup>lt;sup>1</sup> STEDS, Radical Waters (Pty)Ltd., P.O Box 10844 Vorna Valley, 1686, South Africa



Both these solutions remain in the metastable state for approximately 48 hours before the solution returns to the stable state, becoming inactive once more.

ECA is produced from tap water with the addition of saline solution in a special unit (Fig 5) that houses a unique flow-through electrolytic module (FEM).






The FEM contains the anode, made from titanium and coated with ruthenium –oxide, iridium and platinum, and the cathode, made from titanium coated with pyro-carbon and glass-carbon. A diaphragm consisting of ultra-filtration, electro-catalytic ceramics on a bed of zirconium, yttrium, aluminium and niobium-oxides, separates the anode and cathode.

A pilot study was conducted in which the same anaerobic bacteria used in this study were streaked out onto blood agar plates. ECA solutions with different pH values absorbed into blotting paper were placed onto the agar plates and incubated anaerobically at 37°C for 7 days. In a similar manner to antibiotic-sensitivity testing, successful inhibition of growth of all four bacteria for some distance around the ECA

solutions was demonstrated (Fig 6 & 7 – Prevotella intermedia).

Fig 6: Pilot study to demonstrate anti- Fig 7: Higher magnification of blood Microbial properties of ECA

agar plate seen in Fig 6





The research of Hata et al., 1996, reported on the effectiveness of so-called oxidative potential water (OPW) as a root canal irrigant. The authors, from Osaka Dental University, Japan, reported favorably on the cleaning ability of OPW to clean debris from canal walls. It is not exactly clear what OPW is, but there is subjective reason to believe that it may be similar to ECA and being produced by similar technology<sup>2</sup>.

<sup>&</sup>lt;sup>2</sup> Super Mini-Water, Janix Inc., Atsugi, Japan



### **CHAPTER 3**

### Purpose

The purpose of this *in vitro* study was to evaluate the antimicrobial effectiveness of ECA on a selected group of anaerobic bacteria in root canals of extracted human teeth.

### **CHAPTER 4**

### **Materials and Methods**

Sixty caries-free, single rooted, adult, maxillary, anterior human teeth were used for the study (Fig 8).

Fig 8: Teeth used for this study



Teeth were collected from the extraction clinic of the University of Pretoria and rinsed under cold running tap water and stored in distilled water at 4°C for 1 week.

Radiographs were taken of each tooth to identify and eliminate teeth with aberrant canal forms, caries, internal or external resorption, calcifications, multiple canals, or any other condition, which may influence the experimental procedure.



### 4.1 Preparation of teeth

Conventional access cavities were prepared using diamond burs in a turbine driven handpiece and long shank round burs in an air-motor driven handpiece. Canals were explored using number 08 X-files<sup>3</sup>, confirming apical patency and establishing working length visually, for each individual tooth. The working length of each tooth was noted at the time of preparation. The coronal two thirds of the canal was prepared using a nickel-titanium instrument in a rotary handpiece<sup>4</sup>.

After each size file, and intermittently during the preparation procedure, canals were irrigated with a 3.5% solution of NaOCl<sup>5</sup>, delivered from an ultrasonic unit<sup>6</sup>.

The apical third of each canal was prepared to a minimum ISO size 30, using a series of hand-files. Once again, after each size file, the canals were thoroughly irrigated as described above. Following apical preparation, the canals were smoothly tapered using a step-back procedure with X-files and rotary nickel titanium files, again with copious irrigation (NaOCl).

The apices of the teeth were sealed with a resin modified glass-ionomer cement<sup>7</sup>. The remaining root surface was coated with a double layer of nail varnish, to isolate the internal root canal environment (Behrend, Cutler and Gutmann, 1996). Each canal was finally flushed with NaOCl, dried with a paper point and filled with calcium hydroxide<sup>8</sup> (Ca(OH)<sub>2</sub>) powder (Turkun and Cengiz, 1997), mixed with sterile water to a creamy mixture and introduced into the canal with a spiral filler (Fig 9)

<sup>&</sup>lt;sup>3</sup> Hyflex X-File, Hygenic Corporation, Akron, Ohio, USA

<sup>&</sup>lt;sup>4</sup> ProFile Series 29, Maillefer, Ballaigues, Switzerland

<sup>&</sup>lt;sup>5</sup> Jik, Reckitt Colman, Elandsfontein, South Africa

<sup>&</sup>lt;sup>6</sup> Piesoelectric, Siemens, Germany

<sup>&</sup>lt;sup>7</sup> Advance, Caulk/Dentsply, Milford, USA

<sup>&</sup>lt;sup>8</sup> Calcium hydroxide, Merck, D-6100, Darmstadt, Germany



Fig 9: Spiral filler introduction of Ca(OH)<sub>2</sub> (note the sealed external root surface)



A cotton wool pellet was placed into the pulp chamber and the access cavity sealed with a zinc oxide eugenol material<sup>9</sup>. After a 7-day storage period in de-ionized water (37°C), access to the canals was once more established by removing the zinc oxide eugenol material and cotton wool pellets.

The aim of the preparation, irrigation, filling (Ca(OH)<sub>2</sub>) and storage procedure was to produce tapered, sterile root canals. Sterility was maintained for the duration of the experiment, by working in a laminar flow cabinet in a positive sterile airflow laboratory, and using sterilized instrumentation, gloves and masks.

### 4.2 Preparation for first microbiological analysis

The teeth were placed in sterile Ringer's solution (22°C) for five 24-hour cycles, Ringer's solution being replaced with a fresh solution after each cycle. This was done to remove NaOCl and Ca(OH)<sub>2</sub> from the canals and dentinal tubules.

25 microlitre (µl) sterile Ringer's solution was injected, using a sterile pipette, into each canal and aspirated (Fegan and Steiman, 1995)(Fig 10).

<sup>&</sup>lt;sup>9</sup> IRM, Caulk/Dentsply, Milford, USA





This procedure was repeated three times, serving as a rinsing procedure. The aspirated Ringer's solution was placed onto a blood agar plate (Fig 11) and spread,

using a sterile glass rod (Harrigan and McCance, 1970).

samples from root canal







The plates were incubated anaerobically (Anerocult<sup>10</sup>) at 37°C for seven days. All cultures were negative at seven days, confirming absence of anaerobic bacteria within the root canals. Bacterial standardization (Fig 12) was done by using a MacFarland Standard-1 suspension and spectrophotometry<sup>11</sup> at 540 nm.

Anerocult A®, Merck GaA, 64271 Darmstadt, Germany
 GBC -uv/vis 916. Spectrophotometer (Wirsam), GBC Scientific Equipment, Melbourne, Australia

Calibrated suspensions of 2 obligate 12,13 and 2 facultative 14,15 anaerobic bacteria (Fig

13), were prepared and inoculated into 500 ml brain-heart-infusion (BHI) to form a broth (Kwang-ShikBae *et al.*, 1997; Siqueira *et al.*, 1997; Siqueira and de Uzeda, 1996).

# Fig 12; Standardized suspensions of bacteria (Mac Farland-1)



Fig 13: Anaerobic bacteria imported from – ATCC(USA)



<sup>&</sup>lt;sup>12</sup> Prevotella intermedia

ATCC No.25611 ATCC No.33277

ATCC No. 49474

ATCC No.33384

<sup>13</sup> Porphyromonas gingivalis

<sup>&</sup>lt;sup>14</sup> Enterococcus faecalis

<sup>15</sup> Actinobacillus actinomycetemcomitans



Dilution of bacteria (for calibration) was as follows:

 $P \ gingivalis = 133 \times 10^3 \quad P \ intermedia = 5 \times 10^3$ 

E faecalis =  $420 \times 10^3$  A actinomycetemcomitans =  $2 \times 10^3$ 

The morphological appearance of the individual organism and the principles of opacity and absorbency related to spectrophotometry, determine that a smaller bacteria will be inoculated in higher numbers to achieve the same spectrophotometric reading, as for a larger microorganism.

The teeth were placed into a glass beaker containing 500ml BHI broth (inoculated with bacteria) and incubated anaerobically at 37°C, in a shake-incubator, for 7 days. Anaerobic incubation was accomplished by placing the glass beaker into a sealed anaerobic flask containing an Anerocult insert and the flask placed into an incubator. After 7 days the teeth were removed from the broth in an upright position retaining the suspension within the canals (Fig 14).

Fig 14: Tooth removed from BHI broth and held upright (note BHI broth in beaker)





The suspension was aspirated from each canal and 25 µl spread onto marked blood agar plates. The plates were incubated anaerobically at 37°C for seven days and the presence of organisms inside the canal of each tooth was confirmed by formation of colony forming units (exceeding 1000 per plate). Immediately following aspiration, each tooth was placed into a sterile, marked container that corresponded to the marked agar plate for that tooth.

### 4.3 Irrigation of teeth

The teeth were randomly divided into four groups (n=15). All teeth were irrigated for 5 minutes, with the four different irrigation solutions (Table 1:) being delivered from an ultrasonic unit:

**Table 1: Irrigation solutions and Groups** 

Group A	Distilled water (Control group)
(n=15)	(obtained from laboratory)
Group B	Sodium hypochlorite
(n=15)	(Jik – 3.5% concentration) *
Group C	ECA at pH 7.0
(n=15)	(Anolyte) **
Group D	ECA at pH 9.0
(n=15)	(Anolyte) **

<sup>\*</sup> Jik is manufactured at this concentration and is within the range of NaOCl concentrations shown to be antimicrobial (Baumgartner and Quenin, 1992).

### 4.4 Second microbiological analysis

Directly following irrigation, 25µl sterile BHI solution was injected into each canal and aspirated immediately. This procedure was repeated three times (rinsing the canal). The third and final 25µl aspirated BHI solution was inoculated into 5.0ml BHI broth and measured by spectrophotometry.

<sup>\*\*</sup> The solutions for groups C and D were produced by a STEDS unit<sup>16</sup> set at a current of 5.0 Ampere.

<sup>&</sup>lt;sup>16</sup> STEDS, Radical Waters, Vorna Valley, South Africa



In the same manner, a further 25µl BHI solution was aspirated and placed onto a blood agar plate and incubated anaerobically at 37°C for 7 days.

The access cavities were sealed using Fermit<sup>17</sup>, the teeth placed in distilled water and incubated anaerobically at 37°C for 7 days..Fermit was used to facilitate easy removal of the restoration. Canals were accessed after 7 days and spectrophotometric analysis and plate counts for colony forming units repeated.

### 4.5 Statistical analysis

Results of the standard plate counts and spectrophotometry were analyzed using analysis of variance (ANOVA) and Dunnett's one-tailed t-test to detect differences between the four groups and in particular to test for significant differences between the control group (Group A) and the other three groups (Groups B,C and D)

Multiple comparison analysis, using Scheffe, Duncan and Tukey tests, was done to pinpoint differences between groups.

A paired t-test was performed on spectrophotometric readings for each group, to determine whether significant growth of bacteria took place in the period directly following irrigation and the seventh day post-irrigation.

### 4.6 Scanning electron microscopy

Nine teeth were removed at various stages of the experiment (Table 2) and fixed in 2.5% gluteraldehyde and phosphate-buffered solution (0.1M) (Glauert, 1981). Teeth were sectioned longitudinally, using a low-speed (200 rpm), water cooled, diamond saw<sup>18</sup>. All specimens were mounted on metal plates with conductive adhesive paste, critical point dried and sputter coated with an 80 micron thick layer of gold.

Scanning electron micrographs<sup>19</sup> were taken at magnifications 100, 1500 and 6000 of nine teeth (Table 2).

Fermit, Vivadent Ets., FL-9494 Schaan, Lichtenstein
 Isomet, Buehler, Micro-Met Scientific, P.O.Box 84386 Greenside, South Africa
 JEOL 5800, Tokyo, Japan



Table 2: Stages at which teeth were removed from study for SEM examination

Directly from inoculated BHI broth – no irrigation
Following irrigation with distilled water and no incubation
Following irrigation with distilled water and incubation for 7 days
Following irrigation with NaOCl and no incubation
Following irrigation with NaOCl and incubation for 7 days
Following irrigation with ECA pH 7.0 and no incubation.
Following irrigation with ECA pH 7.0 and incubation for 7 days
Following irrigation with ECA pH 9.0 and no incubation
Following irrigation with ECA pH 9.0 and incubation for 7 days

**Note:** Five (of the above nine) teeth were removed from the study before the first microbiological analysis and resulted in sample numbers (n=13), (n=14), (n=14) and (n=14) for the four groups (A, B, C & D) respectively. The remaining four teeth examined by SEM were removed from the study, subsequent to the second microbiological analysis



### **CHAPTER 5**

### **Results**

Results of the standard plate counts (colony forming units) and spectrophotometric analysis are given in Tables 3 and 4.

Table 3: Colony forming units and spectrophotometric values measured immediately after irrigation of root canals with four different solutions

	No. of Teeth	, ,		Spectrophotometric values	
	(n)	<b>Mean</b> (4x10 <sup>1</sup> )/ml	SD	Mean	SD
Group A (Control)	13	>1000*	0	0.6585	0.0959
Group B (NaOCl)	14	0	0	0.0442	0.0134
Group C (ECA pH 7.0)	14	693	253	0.5785	0.1914
Group D (ECA pH 9.0)	14	525	418**	0.6251	0.1377

<sup>\*</sup> A standardized microbiological counting technique was applied (Harrigan and McCance, 1970). Plates were divided into quadrants and a count exceeding 250 for a quadrant was recorded as exceeding 1000 for that plate.

All of the canals in Group A showed more than 1000 colony forming units, indicating the inability of distilled water to destroy bacteria. The canals in Group B however, showed no colony formation and confirm the bactericidal action of NaOCl. Groups C and D show less numbers of colony formation as well as reduced spectrophotometric values compared to the control group (Group A). This reduction however, does not approach the negative (zero) values as recorded for Group B and in fact are closer to those values obtained for the control group (Group A)

<sup>\*\*</sup> A high SD (standard deviation) value was recorded for this group, indicative of a high variation within the group (CV=79.6).



Table 4: Colony forming units and spectrophotometric values measured seven days after irrigation of root canals with four different solutions

	No. of Teeth	Colony form	rming units Spectrophotometric values		otometric
	(n)	<b>Mean</b> (4x10 <sup>1</sup> )/ml	SD	Mean	SD
Group A (Control)	13	>1000*	0	0.8458	0.1138
Group B (NaOCl)	14	114**	291	0.0654	0.0138
Group C (ECA pH 7.0)	14	814	183	0.7524	0.1670
Group D (ECA pH 9.0)	14	721	286	0.8082	0.2467

A standardized microbiological counting technique was applied (Harrigan and McCance, 1970). Plates were divided into four sections and a count exceeding 250 for a quadrant was recorded as exceeding 1000 for that plate.

All four groups show an increase in colony forming units and spectrophotometric values when compared to the results of the immediate post-irrigation analysis (Table 3). Values for Groups C and D once again approximate the values for Group A, and the values for Group B are almost negligible in comparison. If we discard the values of the contaminated plates in the colony forming unit count for Group B, a zero count is again recorded.

<sup>\*\*</sup> Two of the fourteen plates had values of 800 (see Appendix II, pg.29, specimens 20 and 28), resulting in a mean of 114 for the entire group. These colonies could be clearly identified as contaminants.



Table 5: Increase in colony forming units and spectrophotometric values for the period directly following irrigation to the seventh day post-irrigation (Table 4 minus Table 3)

	No. of	Mean ***		
	Teeth (n)	Colony forming units (4x10 <sup>1</sup> )/ml	Spectrophotometric values	
Group A (Control)	13	*	0.1872	
Group B (NaOCl)	14	114**	0.0211	
Group C (ECA Ph 7.0)	14	121	0.1738	
Group D (ECA pH 9.0)	14	196	0.1831	

<sup>\*</sup> This value could not be determined, due to the non-defined values exceeding 1000.

Once more, the values in Groups A, C and D are very close whilst the spectrophotometric value for Group B is at least 10<sup>1</sup> less than that of the other groups and the number of colonies formed is again zero (disregarding contaminants)

<sup>\*\*</sup> Two of the fourteen plates had values of 800, resulting in a mean of 114 for the entire group. These colonies could be clearly identified as contaminants.

<sup>\*\*\*</sup> Standard deviation values are not given, as these were not statistically representative and have no meaningful interpretation.



Table 6: Dunnett's test results on spectrophotometric values, comparing each of Groups B, C and D to Group A (Control)

Comparison	Directly following irrigation	Seven days post-irrigation
Group B to Control	***	***
Group C to Control	-	-
Group D to Control	-	-

<sup>\*\*\*</sup> denotes significant difference (p  $\leq 0.05$ )

Spectrophotometric values for Group B differed significantly from Group A ( $p \le 0.05$ ). The close approximation of values for Groups A, C and D was also demonstrated by this test (p > 0.05), (confirming what was seen in Tables 3, 4 and 5).

Table 7: Results of paired t-test on spectrophotometric values for bacterial growth over the seven day period after irrigation with four different solutions

Group	Mean	p-value	
Group A	0.1872	0.0001*	
(Control)			
Group B	0.0211**	0.0001*	············
(NaOCl)			
Group C	0.1738	0.0027*	
(ECA Ph 7.0)			
Group D	0.1831	0.0015*	
(ECA pH 9.0)			

<sup>\*</sup> Significant microbial growth was found in all four groups for the seven day period between the measuring of microbial counts (p ≤0.05)

<sup>\*\*</sup> Bacterial growth within the canals irrigated with NaOCl was at least one thousandth less than that in any of the other three groups.



### **SEM Examination**

Representative micrographs of the surfaces of the root canals are shown in Figures 15 to 22

Fig 15: No Irrigation (x 6000)




Fig 15: A layer of bacteria was seen on the entire surface of the canal in the control group (Group A), before irrigation was performed. (Original magnification x 6000; bar indicates 1 micron)

Fig 16: Distilled water (x 6000)




Fig 16: Surface of a root canal of a tooth in Group A following irrigation with distilled water. Abundant debris, presumably covering bacteria was seen. (Original magnification x 6000; bar indicates 1 micron)

Fig 17: NaOCl (x 1500)




Fig 17: Surface of root canal of a tooth in Group B irrigated with NaOCl. Note the presence of the smear layer partially occluding dentinal tubuli. (Original magnification x 1500; bar indicates 10 micron)

Fig 18: NaOCl (x 6000)




Fig 18: Higher magnification of area shown in Fig.17 (Original magnification x 6000; bar indicates 1 micron)

Fig 19: ECA pH 7.0 (x 1500)




Fig 19: Surface of root canal of a tooth in Group C irrigated with ECA, pH 7.0. Lack of debris is evident. (Original magnification x 1500; bar indicates 10 micron)

Fig 20: ECA pH 7.0 (x 6000)




Fig 20: Higher magnification of area shown in Fig.19, showing thinner smear layer than that seen in Fig.18 (NaOCl). (Original magnification x 6000; bar indicates 1 micron)

Fig 21: ECA pH 9.0 (x 1500)




Fig 21: Surface of root canal of a tooth in Group D irrigated with ECA, pH 9.0 showing lack of debris or bacteria. (Original magnification x 1500; bar indicates 10 micron)

Fig 22: ECA pH 9.0 (x 6000)




Fig 22: Higher magnification of root canal of a tooth in Group D irrigated with ECA, pH 9.0 (different tooth and root canal to that shown in Fig. 21 above), showing removal of smear layer. (Original magnification x 6000; bar indicates 1 micron)



### **CHAPTER 6**

### **Discussion**

The selection of the four bacteria used in this study was based on the fact that they are all widely implicated pathogens in diseases of the pulp and periapical tissues. In addition, these bacteria have also been used, either as individual strains or in combination, in many *in vitro* studies of this nature (Safavi *et al*, 1990; Yesilsoy *et al*, 1995; Siqueira and de Uzeda, 1996; Siren *et al*, 1997; Siqueira *et al*, 1997; Kwang-Shik *et al*, 1997; Barbosa *et al*, 1997; Siqueira and de Uzeda, M, 1997; Dougherty *et al*, 1998; Molander *et al*, 1998).

Taxonomic identification was not carried out for growth on the culture medium and thus any and all remaining bacteria were included in the colony forming unit count. Because the study was conducted on a limited number of pure cultured (ATCC, USA) anaerobic bacteria and under sterile conditions, there was no intention or need to do taxonomic identification.

Contamination was apparent in only two of the total number of 240 plates cultured and this was easily identified (green, mucous patches) and did not affect our findings. Of interest, is the fact that this contamination was found in the group irrigated with NaOCl.

All of the cultures for the control group had bacterial counts exceeding 1000 colony forming units, representative of active bacterial growth. The fact that these bacteria successfully withstood vigorous irrigation (ultrasonics) with distilled water and to a lesser extent, irrigation with ECA, proves that the four selected species of bacteria are a useful model for a study of this kind.

The fact that colony forming units consistently exceeded 1000 (undefined) in the control group (Group A), was reason to consider a larger dilution factor (i.e.  $> 4x10^1$ ). Should defined counts have been possible, this would have made statistical evaluation even more accurate.



However, the probability of recording 'false-negative' (decreased sensitivity for the test) counts for Groups C and D (see Appendix I-specimens 48, 54, 55, 56 and 59) would have been more detrimental to the results and the dilution factor thus used was  $4x10^{1}$  (0.25µl)

Under the conditions of this study, ECA failed to destroy all of the bacteria within the root canals of the teeth studied. On the other hand NaOCl reduced the bacterial counts to less than 1 colony forming unit /  $25\mu l$  in all of the specimens. This was the case both directly following irrigation and very importantly, also seven days post-irrigation. Table 7 shows that although there was an increase in the number of bacteria present in the root canals over the seven day incubation period, Group B (NaOCl) had by far the least (Mean = 0.0211). This confirms the total antimicrobial action of NaOCl.

The control group consistently (15 from 15 counts at each analysis) produced bacterial counts greater than 1000 colony forming units/ 25µl, confirming the accuracy and usefulness of this model for similar studies. The results of the spectrophotometric analysis correlated well with those for the colony forming units (as seen in Tables 3, 4 and 5), further confirming the relevance of this model. Although subjective, visual examination by SEM supported these findings by showing many bacteria in the control group and no bacteria in the group irrigated with NaOCl. Varying numbers of bacteria were seen for the groups irrigated with ECA (Fig.15-22).

The large variation (SD) recorded for colony forming units for Group D (Table 3) could not be accounted for. Thorough investigation as to the possible reason for this occurrence revealed no explanation. The very low (groups A and B) and low (Group C) variation, of all three other groups, seem to indicate that the experimental model and techniques used were accurate. This, together with the fact that values for spectrophotometry indicate an acceptable level of variation (2SD) in all four groups and furthermore, the fact that this high variation was not seen in the measurements taken seven days post-irrigation, give further reason to accept this model as accurate.



Although fewer bacteria were measured and seen (SEM) in the groups irrigated with ECA (compared to Group A, control) both directly following irrigation and seven days post-irrigation, ECA was not as effective in destroying bacteria as was NaOCl. Various reasons for this apparent failure may be given. Firstly, there may be an inherent problem with the technological production of ECA, however, in the light of the successful pilot studies conducted prior to commencement of this study, this is doubtful. Secondly, the results of earlier work on ECA (Marais, 1998) suggest that eradication of bacteria within the root canal system may be achieved with the use of ECA, another inconsistency with the apparent failure of ECA in this study.

A more likely explanation may be one of the following:

- 1. At the time of preparation of the ECA solutions used in this study, a technical failure of the unit may have resulted in an altered ECA solution. This would indicate the technique sensitivity in the production of ECA for use as an endodontic solution and indicates the need for quality control of the technology.
- 2. The specific ECA parameters (pH, current, temperature, exposure time of canal to irrigant, flow rate) used in this study may need to be varied and will be studied in future research. As can be expected in any system with so many variables, and given the metastable state of the solution, one can expect to have initial difficulties in identifying the exact parameters required for a successful antimicrobial solution for use in the root canal system.
- 3. In the study by Marais (1998), the irrigation procedure with the anolyte was preceded by a similar rinse with catholyte. In the present study, only anolyte was used. The use of anolyte and catholyte in combination, should be investigated further.



### **CHAPTER 7**

### **Conclusion**

- 1. In accordance with many other studies, NaOCl (3,5%) was shown to be completely bactericidal within the root canal system.
- 2. The use of ECA caused a reduction in the number of anaerobic bacteria within the root canal system, but this was not statistically significant (p>0.05).
- 3. The model used for this study is proposed for future research into antimicrobial effectiveness against anaerobic bacteria in the root canal system.

"It is not enough that you should understand about applied science in order that your work may increase man's blessings. Concern for Man himself and his fate must always form the chief interest of all technical endeavors.... In order that the creations of our mind shall be a blessing and not a curse to mankind – never forget this in the midst of your diagrams and equations." Albert Einstein



### REFERENCES

American Association of Endodontics recertification document, (1989)

Barbosa CAM, Goncalves RB, Siqueira JF, de Uzeda M (1997) Evaluation of the antibacterial activities of calcium hydroxide, chlorhexidine, and camphorated paramonochlorophenol as intracanal medicament. A clinical and laboratory study. *Journal Of Endodontics* **23**, 297-299.

Baumgartner JC, Cuenin PR (1992) Efficacy of several concentrations of sodium hypochlorite for root canal irrigation. *Journal Of Endodontics* **18**, 605-612.

Behrend GD, Cutler CW, Gutmann JL (1996) An *in vitro* study of smear layer removal and microbial leakage along root-canal fillings. *International Endodontic Journal* **29**, 99-107.

Brown DC, Moore BK, Brown CE Jr., Newton CW (1995) An *in vitro* study of apical extrusion of sodium hypochlorite during endodontic canal preparation. *Journal of Endodontics* **21**, 587-591.

Caliskan MK, Turkun M, Alper S (1994) Allergy to sodium hypochlorite during root canal therapy: a case report. *International Endodontic Journal* **27,** 163-167.

Cameron JA (1995) Factors affecting the clinical efficiency of ultrasonic endodontics: a scanning electron microscopy study. *International Endodontic Journal* **28**, 47-53.

Cheung GS, Stock CJR (1993) *In vitro* cleaning ability of root canal irrigants with and without endosonics. *International Endodontic Journal* **26**, 334-343.

Cymbler DM, Ardakani P (1994) Sodium hypochlorite injection into periapical tissues. *Dental Update* 21, 345-346.



Dougherty WJ, Bae KS, Watkins MS, Baumgartner JC (1998) Black-pigmented bacteria in coronal and apical segments of infected root canals. *Journal of Endodontics* **24**, 356-358.

Fegan SE, Steiman HR (1995) Comparative evaluation of the antibacterial effects of intracanal Nd-YAG laser irradiation: an *in vitro* study. *Journal of Endodontics* **21**, 415-417.

Georgopoulou M, Kontakiotis E, Nakou M (1994) Evaluation of the antimicrobial effectiveness of citric acid and sodium hypochlorite on the anaerobic flora of the infected root canal. *International Endodontic Journal* **27**, 139-143.

Glauert AM (1981) Fixation, Dehydration and Embedding of Biological Specimens. In: *Practical Methods In Electron Microscopy*: 4<sup>th</sup> Edition, North Holland Publishing Company, Amsterdam.

Gulabivala K, Stock CJR (1995) Preparation of the root canal. In: *Color Atlas and Text of Endodontics*. 2<sup>nd</sup> Edition, Mosby-Wolfe, 97-144.

Harrigan WF, McCance ME (1970) Laboratory Methods in Microbiology. Part1:34-35. London Academic Press.

Hata G-I, Uemura M, Weine F, Toda T (1996) Removal of smear layer in the root canal using oxidative potential water. *Journal of Endodontics* **22**, 643-645.

Jeansonne MJ, White RR (1994) A comparison of 2% chlorhexidine gluconate and 5,25% sodium hypochlorite as antimicrobial endodontic irrigants. *Journal of Endodontics* **20**, 276-278.

Jimenez-Rubio A, Segura JJ, Llamas R, Jimenez-Planas A, Guerrero JM, Calvo JR (1997) *In vitro* study of the effect of sodium hypochlorite and gluteraldehyde on substrate adherence capacity of macrophages. *Journal of Endodontics* **23,**562-564



Kwang-Shik Bae, Baugartner JC, Shearer TR, David LL (1997) Occurrence of *Prevotella nigrescens* and *Prevotella intermedia* in infections of endodontic origin. *Journal of Endodontics* 23, 620-623.

Marais JT (1998) Cleaning efficacy of a new root canal irrigation material Journal of Dental Research 77, Abstract 300, 669.

Molander A, Reit C, Dahlen G, Kvist, T (1998) Microbiological status of root-filled teeth with apical periodontitis. *International Endodontic Journal* **31**, 1-7.

Safavi K, Spangberg L, Langeland K (1990) Root canal dentinal tubule disinfection. Journal of Endodontics 16, 207-210.

Siqueira JF, Machado AG, Silveira RM, Lopes HP, de Uzeda, M (1997) Evaluation of the effectiveness of sodium hypochlorite used with three irrigation methods in the elimination of *Enterococcus faecalis* from the root canal, in vitro. *International Endodontic Journal* 30, 279-282.

Siqueira JF, de Uzeda, M (1996) Disinfection by calcium hydroxide pastes of dentinal tubules infected with two obligate and one facultative anaerobic bacteria. *Journal of Endodontics* **22**, 674-676.

Siqueira, JF, de Uzeda M (1997) Intracanal medicaments: Evaluation of the antibacterial effects of chlorhexidine, metronidazole, and calcium hydroxide associated with three vehicles. *Journal of Endodontics* **23**, 167-169.

Siren EK, Haapasalo MPP, Ranta K, Salmi P, Kerusuo E.N.J (1997) Microbiological findings and clinical treatment procedures in endodontic cases selected for microbiological investigation. *International Endodontic Journal* **30**, 91-95.

Turkun M, Cengiz, T (1997) The effects of sodium hypochlorite and calcium hydroxide on tissue dissolution and root canal cleanliness. *International Endodontic Journal* **30**, 335-342.



Walton RE, Rivera, EM (1996) Cleaning and shaping. In: *Principles and Practice of Endodontics*. 2<sup>nd</sup> Edition WB Saunders, Philadelphia 201-233.

Yesilsoy C, Whitaker E, Cleveland D, Phillips E, Trope M (1995) Anti-microbial and toxic effects of established and potential root canal irrigants. *Journal of Endodontics* **21**, 513-515.



# **ADDENDA**

### Appendix A

Colony forming units and spectrophotometric analysis measured directly after Irrigation

Number units (4x10¹)/ml         analysis units (4x10¹)/ml         Number (4x10¹)/ml         forming units (4x10¹)/ml         analysis units (4x10¹)/ml           1         >1000         0,5504         31         XXXXXXX         XXXXXXX           2         >1000         0,5851         32         >1000         0,5690           3         >1000         0,6167         34         >800         0,5185           5         >1000         0,6193         35         >500         0,8329           6         XXXXXX         XXXXXXX         36         >1000         0,5773           7         XXXXXXX         XXXXXXXX         37         >500         0,5639           8         >1000         0,6080         39         >500         0,6995           9         >1000         0,6080         39         >500         0,4952           10         >1000         0,7985         40         >1000         0,5559           11         >1000         0,6096         41         >500         0,5744           12         >1000         0,7083         42         >500         0,5549           13         >1000         0,8843         43         >500         0,7138	irrigation					
units (4x10¹)/ml         units (4x10¹)/ml           1         >1000         0,5504         31         XXXXXX         XXXXXXX           2         >1000         0,5851         32         >1000         0,5696           3         >1000         0,6386         33         >1000         0,6045           4         >1000         0,6167         34         >800         0,5185           5         >1000         0,6193         35         >500         0,8329           6         XXXXXX         XXXXXXX         36         >1000         0,5773           7         XXXXXX         XXXXXXX         37         >500         0,5639           8         >1000         0,7340         38         >500         0,6995           9         >1000         0,6080         39         >500         0,4952           10         >1000         0,7985         40         >1000         0,5559           11         >1000         0,6086         41         >500         0,5549           12         >1000         0,7083         42         >500         0,5549           13         >1000         0,5860         44         >400         0	Specimen	Colony	Spectrometric	Specimen	Colony	Spectrometric
(4x10¹)/ml         (4x10¹)/ml           1         >1000         0,5504         31         XXXXXX         XXXXXXX           2         >1000         0,5851         32         >1000         0,5690           3         >1000         0,6386         33         >1000         0,6045           4         >1000         0,6167         34         >800         0,5185           5         >1000         0,6193         35         >500         0,8329           6         XXXXXX         XXXXXXX         36         >1000         0,5773           7         XXXXXX         XXXXXXX         37         >500         0,5639           8         >1000         0,7340         38         >500         0,6995           9         >1000         0,6080         39         >500         0,4952           10         >1000         0,7985         40         >1000         0,5559           11         >1000         0,7983         42         >500         0,5549           13         >1000         0,7883         42         >500         0,5549           13         >1000         0,5960         44         >400         0,0473	Number	forming	analysis	Number	forming	analysis
1         >1000         0,5504         31         XXXXXX         XXXXXXX           2         >1000         0,5851         32         >1000         0,5690           3         >1000         0,6386         33         >1000         0,6045           4         >1000         0,6167         34         >800         0,5185           5         >1000         0,6193         35         >500         0,8329           6         XXXXXX         XXXXXXXX         36         >1000         0,5773           7         XXXXXX         XXXXXXXX         37         >500         0,5639           8         >1000         0,7340         38         >500         0,6995           9         >1000         0,6080         39         >500         0,6995           9         >1000         0,6080         39         >500         0,4952           10         >1000         0,6086         39         >500         0,4952           11         >1000         0,6096         41         >500         0,5579           11         >1000         0,7083         42         >500         0,5549           13         >1000         0,						
2         >1000         0,5851         32         >1000         0,5690           3         >1000         0,6386         33         >1000         0,6045           4         >1000         0,6167         34         >800         0,5185           5         >1000         0,6193         35         >500         0,8329           6         XXXXXXX         XXXXXXX         36         >1000         0,573           7         XXXXXXX         XXXXXXX         37         >500         0,5639           8         >1000         0,7340         38         >500         0,6995           9         >1000         0,6080         39         >500         0,6995           9         >1000         0,6080         39         >500         0,4952           10         >1000         0,7985         40         >1000         0,5559           11         >1000         0,6096         41         >500         0,5074           12         >1000         0,7083         42         >500         0,5549           13         >1000         0,8843         43         >500         0,7138           14         >1000         0,59		$(4x10^{1})/ml$			$(4x10^{1})/ml$	
3         >1000         0,6386         33         >1000         0,6045           4         >1000         0,6167         34         >800         0,5185           5         >1000         0,6193         35         >500         0,8329           6         XXXXXX         XXXXXXX         36         >1000         0,5773           7         XXXXXXX         XXXXXXXX         37         >500         0,5639           8         >1000         0,7340         38         >500         0,6995           9         >1000         0,6080         39         >500         0,4952           10         >1000         0,7985         40         >1000         0,5559           11         >1000         0,7083         42         >500         0,5549           13         >1000         0,7083         42         >500         0,5549           13         >1000         0,8843         43         >500         0,7138           14         >1000         0,5960         44         >400         0,0473           15         >1000         0,6129         45         >1000         0,7656           18         0         0,037		>1000	0,5504		XXXXXX	XXXXXXX
4         >1000         0,6167         34         >800         0,5185           5         >1000         0,6193         35         >500         0,8329           6         XXXXXX         XXXXXXX         36         >1000         0,5773           7         XXXXXX         XXXXXXX         37         >500         0,5639           8         >1000         0,7340         38         >500         0,6995           9         >1000         0,6080         39         >500         0,4952           10         >1000         0,6086         39         >500         0,4952           10         >1000         0,7985         40         >1000         0,5559           11         >1000         0,6096         41         >500         0,5549           13         >1000         0,7083         42         >500         0,5549           13         >1000         0,8843         43         >500         0,7138           14         >1000         0,5960         44         >400         0,0473           15         >1000         0,6129         45         >1000         0,8599           16         XXXXXX         XX			0,5851	32	>1000	0,5690
5         >1000         0,6193         35         >500         0,8329           6         XXXXXX         XXXXXXX         36         >1000         0,5773           7         XXXXXX         XXXXXXX         37         >500         0,5639           8         >1000         0,7340         38         >500         0,6995           9         >1000         0,6080         39         >500         0,4952           10         >1000         0,7985         40         >1000         0,5559           11         >1000         0,6096         41         >500         0,5074           12         >1000         0,7083         42         >500         0,5549           13         >1000         0,8843         43         >500         0,7138           14         >1000         0,5960         44         >400         0,0473           15         >1000         0,6129         45         >1000         0,8599           16         XXXXXX         XXXXXXX         XXXXXXX         XXXXXXX         XXXXXXX           17         0         0,0429         47         >1000         0,7656           18         0 <t< td=""><td></td><td>&gt;1000</td><td>0,6386</td><td></td><td>&gt;1000</td><td></td></t<>		>1000	0,6386		>1000	
6         XXXXXX         XXXXXXX         36         >1000         0,5773           7         XXXXXX         XXXXXXX         37         >500         0,5639           8         >1000         0,7340         38         >500         0,6995           9         >1000         0,6080         39         >500         0,4952           10         >1000         0,7985         40         >1000         0,5559           11         >1000         0,6096         41         >500         0,5074           12         >1000         0,7083         42         >500         0,5549           13         >1000         0,8843         43         >500         0,7138           14         >1000         0,5960         44         >400         0,0473           15         >1000         0,6129         45         >1000         0,8599           16         XXXXXXX         XXXXXXX         XXXXXXX         XXXXXXX         XXXXXXX           17         0         0,0429         47         >1000         0,7656           18         0         0,0371         48         46         0,7507           19         0         0,	4	>1000	0,6167	34	>800	0,5185
7         XXXXXX         XXXXXXX         37         >500         0,5639           8         >1000         0,7340         38         >500         0,6995           9         >1000         0,6080         39         >500         0,4952           10         >1000         0,7985         40         >1000         0,5559           11         >1000         0,6096         41         >500         0,5074           12         >1000         0,7083         42         >500         0,5549           13         >1000         0,8843         43         >500         0,7138           14         >1000         0,5960         44         >400         0,0473           15         >1000         0,6129         45         >1000         0,8599           16         XXXXXXX         XXXXXXX         XXXXXXX         XXXXXXX         XXXXXXX           17         0         0,0429         47         >1000         0,7656           18         0         0,0371         48         46         0,7507           19         0         0,0461         49         >300         0,5529           20         0         0,0679 </td <td>5</td> <td>&gt;1000</td> <td>0,6193</td> <td>35</td> <td>&gt;500</td> <td>0,8329</td>	5	>1000	0,6193	35	>500	0,8329
8         >1000         0,7340         38         >500         0,6995           9         >1000         0,6080         39         >500         0,4952           10         >1000         0,7985         40         >1000         0,5559           11         >1000         0,6096         41         >500         0,5074           12         >1000         0,7083         42         >500         0,5549           13         >1000         0,8843         43         >500         0,7138           14         >1000         0,5960         44         >400         0,0473           15         >1000         0,6129         45         >1000         0,8599           16         XXXXXX         XXXXXXX         XXXXXXX         XXXXXXX         XXXXXXX           17         0         0,0429         47         >1000         0,7656           18         0         0,0371         48         46         0,7507           19         0         0,0461         49         >300         0,5529           20         0         0,0679         50         >1000         0,4822           21         0         0,0426	6	XXXXXX	XXXXXXX	36	>1000	0,5773
9         >1000         0,6080         39         >500         0,4952           10         >1000         0,7985         40         >1000         0,5559           11         >1000         0,6096         41         >500         0,5074           12         >1000         0,7083         42         >500         0,5549           13         >1000         0,8843         43         >500         0,7138           14         >1000         0,5960         44         >400         0,0473           15         >1000         0,6129         45         >1000         0,8599           16         XXXXXX         XXXXXXX         XXXXXXX         XXXXXXXX           17         0         0,0429         47         >1000         0,7656           18         0         0,0371         48         46         0,7507           19         0         0,0461         49         >300         0,5529           20         0         0,0679         50         >1000         0,4822           21         0         0,0426         51         >800         0,5688           22         0         0,0458         52		XXXXXX	XXXXXXX	37	>500	0,5639
10         >1000         0,7985         40         >1000         0,5559           11         >1000         0,6096         41         >500         0,5074           12         >1000         0,7083         42         >500         0,5549           13         >1000         0,8843         43         >500         0,7138           14         >1000         0,5960         44         >400         0,0473           15         >1000         0,6129         45         >1000         0,8599           16         XXXXXX         XXXXXXX         XXXXXXX         XXXXXXX         XXXXXXX         XXXXXXX           17         0         0,0429         47         >1000         0,7656           18         0         0,0371         48         46         0,7507           19         0         0,0461         49         >300         0,5529           20         0         0,0679         50         >1000         0,4822           21         0         0,0426         51         >800         0,7222           23         0         0,0396         53         >1000         0,5539           24         0	8	>1000	0,7340	38	>500	0,6995
11         >1000         0,6096         41         >500         0,5074           12         >1000         0,7083         42         >500         0,5549           13         >1000         0,8843         43         >500         0,7138           14         >1000         0,5960         44         >400         0,0473           15         >1000         0,6129         45         >1000         0,8599           16         XXXXXX         XXXXXXX         46         XXXXXX         XXXXXXX           17         0         0,0429         47         >1000         0,7656           18         0         0,0371         48         46         0,7507           19         0         0,0461         49         >300         0,5529           20         0         0,0679         50         >1000         0,4822           21         0         0,0426         51         >800         0,5688           22         0         0,0458         52         >800         0,7222           23         0         0,0396         53         >1000         0,5539           24         0         0,0305         54 <td>9</td> <td>&gt;1000</td> <td>0,6080</td> <td>39</td> <td>&gt;500</td> <td>0,4952</td>	9	>1000	0,6080	39	>500	0,4952
12         >1000         0,7083         42         >500         0,5549           13         >1000         0,8843         43         >500         0,7138           14         >1000         0,5960         44         >400         0,0473           15         >1000         0,6129         45         >1000         0,8599           16         XXXXXX         XXXXXXX         XXXXXXX         XXXXXXX         XXXXXXXX           17         0         0,0429         47         >1000         0,7656           18         0         0,0371         48         46         0,7507           19         0         0,0461         49         >300         0,5529           20         0         0,0679         50         >1000         0,4822           21         0         0,0426         51         >800         0,5688           22         0         0,0458         52         >800         0,7222           23         0         0,0396         53         >1000         0,5539           24         0         0,0305         54         30         0,5621           25         0         0,0110         55 </td <td>10</td> <td>&gt;1000</td> <td>0,7985</td> <td>40</td> <td>&gt;1000</td> <td>0,5559</td>	10	>1000	0,7985	40	>1000	0,5559
13         >1000         0,8843         43         >500         0,7138           14         >1000         0,5960         44         >400         0,0473           15         >1000         0,6129         45         >1000         0,8599           16         XXXXXX         XXXXXXX         46         XXXXXX         XXXXXXX           17         0         0,0429         47         >1000         0,7656           18         0         0,0371         48         46         0,7507           19         0         0,0461         49         >300         0,5529           20         0         0,0679         50         >1000         0,4822           21         0         0,0426         51         >800         0,5688           22         0         0,0458         52         >800         0,7222           23         0         0,0396         53         >1000         0,5539           24         0         0,0305         54         30         0,5621           25         0         0,0110         55         92         0,4924           26         0         0,0561         56         <	11	>1000	0,6096	41	>500	0,5074
14         >1000         0,5960         44         >400         0,0473           15         >1000         0,6129         45         >1000         0,8599           16         XXXXXX         XXXXXXX         46         XXXXXX         XXXXXXX           17         0         0,0429         47         >1000         0,7656           18         0         0,0371         48         46         0,7507           19         0         0,0461         49         >300         0,5529           20         0         0,0679         50         >1000         0,4822           21         0         0,0426         51         >800         0,5688           22         0         0,0458         52         >800         0,7222           23         0         0,0396         53         >1000         0,5539           24         0         0,0305         54         30         0,5621           25         0         0,0110         55         92         0,4924           26         0         0,0561         56         77         0,5840           27         0         0,0471         57         >80	12	>1000	0,7083	42	>500	0,5549
15         >1000         0,6129         45         >1000         0,8599           16         XXXXXX         XXXXXXX         46         XXXXXX         XXXXXXX           17         0         0,0429         47         >1000         0,7656           18         0         0,0371         48         46         0,7507           19         0         0,0461         49         >300         0,5529           20         0         0,0679         50         >1000         0,4822           21         0         0,0426         51         >800         0,5688           22         0         0,0458         52         >800         0,7222           23         0         0,0396         53         >1000         0,5539           24         0         0,0305         54         30         0,5621           25         0         0,0110         55         92         0,4924           26         0         0,0561         56         77         0,5840           27         0         0,0471         57         >800         0,9350           28         0         0,0554         58         >300 <td>13</td> <td>&gt;1000</td> <td>0,8843</td> <td>43</td> <td>&gt;500</td> <td>0,7138</td>	13	>1000	0,8843	43	>500	0,7138
16         XXXXXX         XXXXXXX         46         XXXXXX         XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX	14	>1000	0,5960	44	>400	0,0473
16         XXXXXX         XXXXXXX         46         XXXXXX         XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX	15	>1000	0,6129	45	>1000	0,8599
$\begin{array}{c ccccccccccccccccccccccccccccccccccc$	16	XXXXXX	XXXXXXX	46	XXXXXX	
19         0         0,0461         49         >300         0,5529           20         0         0,0679         50         >1000         0,4822           21         0         0,0426         51         >800         0,5688           22         0         0,0458         52         >800         0,7222           23         0         0,0396         53         >1000         0,5539           24         0         0,0305         54         30         0,5621           25         0         0,0110         55         92         0,4924           26         0         0,0561         56         77         0,5840           27         0         0,0471         57         >800         0,9350           28         0         0,0554         58         >300         0,5249           29         0         0,0413         59         103         0,4898	17	0	0,0429	47	>1000	0,7656
20         0         0,0679         50         >1000         0,4822           21         0         0,0426         51         >800         0,5688           22         0         0,0458         52         >800         0,7222           23         0         0,0396         53         >1000         0,5539           24         0         0,0305         54         30         0,5621           25         0         0,0110         55         92         0,4924           26         0         0,0561         56         77         0,5840           27         0         0,0471         57         >800         0,9350           28         0         0,0554         58         >300         0,5249           29         0         0,0413         59         103         0,4898	18	0	0,0371	48	46	0,7507
21         0         0,0426         51         >800         0,5688           22         0         0,0458         52         >800         0,7222           23         0         0,0396         53         >1000         0,5539           24         0         0,0305         54         30         0,5621           25         0         0,0110         55         92         0,4924           26         0         0,0561         56         77         0,5840           27         0         0,0471         57         >800         0,9350           28         0         0,0554         58         >300         0,5249           29         0         0,0413         59         103         0,4898	19	0	0,0461	49	>300	0,5529
22         0         0,0458         52         >800         0,7222           23         0         0,0396         53         >1000         0,5539           24         0         0,0305         54         30         0,5621           25         0         0,0110         55         92         0,4924           26         0         0,0561         56         77         0,5840           27         0         0,0471         57         >800         0,9350           28         0         0,0554         58         >300         0,5249           29         0         0,0413         59         103         0,4898	20	0	0,0679	50	>1000	0,4822
23         0         0,0396         53         >1000         0,5539           24         0         0,0305         54         30         0,5621           25         0         0,0110         55         92         0,4924           26         0         0,0561         56         77         0,5840           27         0         0,0471         57         >800         0,9350           28         0         0,0554         58         >300         0,5249           29         0         0,0413         59         103         0,4898	21	0	0,0426	51	>800	0,5688
24         0         0,0305         54         30         0,5621           25         0         0,0110         55         92         0,4924           26         0         0,0561         56         77         0,5840           27         0         0,0471         57         >800         0,9350           28         0         0,0554         58         >300         0,5249           29         0         0,0413         59         103         0,4898	22	0	0,0458	52	>800	0,7222
25     0     0,0110     55     92     0,4924       26     0     0,0561     56     77     0,5840       27     0     0,0471     57     >800     0,9350       28     0     0,0554     58     >300     0,5249       29     0     0,0413     59     103     0,4898	23	0	0,0396	53	>1000	0,5539
26     0     0,0561     56     77     0,5840       27     0     0,0471     57     >800     0,9350       28     0     0,0554     58     >300     0,5249       29     0     0,0413     59     103     0,4898	24	0	0,0305	54	30	0,5621
27     0     0,0471     57     >800     0,9350       28     0     0,0554     58     >300     0,5249       29     0     0,0413     59     103     0,4898	25	0	0,0110	55	92	0,4924
28     0     0,0554     58     >300     0,5249       29     0     0,0413     59     103     0,4898	26	0	0,0561	56	77	0,5840
28     0     0,0554     58     >300     0,5249       29     0     0,0413     59     103     0,4898	27	0			>800	
29 0 0,0413 59 103 0,4898	28	0		58		
	29	0				
	30	0	0,0565	60	>1000	0,7675

<sup>&#</sup>x27;XXX' denotes teeth for SEM evaluation



# Appendix B

Colony forming units and spectrophotometric analysis measured seven days after irrigation

after irrigation					
Specimen	Colony	Spectrometric	Specimen	Colony	Spectrometric
Number	forming	analysis	Number	forming	analysis
	units			units	
	$(4x10^1)/ml$			$(4x10^1)/ml$	
1	>1000	0,6980	31	XXXXXX	XXXXXXX
2	>1000	0,6901	32	>1000	0,6687
3	>1000	0,8339	33	>1000	0,7445
4	>1000	0,7865	34	>1000	0,6329
5	>1000	0,8976	35	>800	0,8936
6	XXXXXX	XXXXXXX	36	>1000	0,6715
7	XXXXXX	XXXXXXX	37	>600	0,5998
8	>1000	0,8310	38	>600	0,7002
9	>1000	0,9802	39	>800	0,6653
10	>1000	1,0140	40	>1000	0,7994
11	>1000	0,7834	41	>600	0,7574
12	>1000	0,8511	42	>600	0,6182
13	>1000	0,9971	43	>800	0,8077
14	>1000	0,9354	44	>600	0,7155
15	>1000	0,6981	45	>1000	1,2590
16	XXXXXX	XXXXXXX	46	XXXXXX	XXXXXXX
17	0	0,0852	47	>1000	0,8551
18	0	0,0512	48	>600	0,8147
19	0	0,0653	49	>800	0,6345
20	>800	0,0799	50	>1000	0,8686
21	0	0,0532	51	>800	0,7347
22	0	0,0776	52	>800	1,3772
23	0	0,0597	53	>1000	0,7598
24	0	0,0411	54	>400	0,6735
25	0	0,0451	55	>300	0,5617
26	0	0,0774	56	>300	0,6250
27	0	0,0617	57	>1000	1,2942
28	>800	0,0762	58	>800	0,6440
29	0	0,0683	59	>300	0,5981
30	0	0,0741	60	>1000	0,8747
/					

<sup>&#</sup>x27;XXX' denotes teeth for SEM evaluation



### Appendix C

The anolyte produced at the anode contains  $Cl_2$ , HClO, ClO, ClO, Cl,  $HO_2$ ,  $HO_2$ ,  $O_2$ , HO,  $O_3$ ,  $O_2$ ,  $^3O_2$ ,  $^1O_2$ , O,  $H_3O^+$ , Cl, H,  $H_2O_2$ ,  $Cl_2O$ ,  $ClO_2^-$ , HCl,  $Cl_2O_7$ ,  $S_2O_8^{8^-}$ ,  $C_2O_6^{2^-}$ , HclO,  $H_2SO_4$ , and  $HSO_3Cl$ .

The catholyte produced at the cathode contains  $HO^-$ ,  $H_3O_2^-$ ,  $O_2^-$ ,  $HO_2^-$ ,  $H_2O_2$ ,  $H_2$ ,  $HO_2$ 

(Radical Waters (Pty)Ltd., P.O Box 10844 Vorna Valley, 1686, South Africa)