

**Table S7. Dehydration and wax infiltration of stem samples.** Dehydration solutions were diluted in distilled water.

<b>Dehydration solution</b>	<b>Incubation</b>
25% butanol, 30% ethanol	60 minutes
40% butanol, 30% ethanol	60 minutes
55% butanol, 25% ethanol	60 minutes
70% butanol, 20% ethanol	120 minutes
85% butanol, 15% ethanol	120 minutes
100% butanol	60 minutes
<b>Wax infiltration</b>	
50% wax, 50% butanol	Overnight, 60 °C
100% wax	Until butanol evaporates, 60 °C